

Updated: Jul 3, 2024

Adapted from the Click-iT EdU Imaging Kit protocol ([C10340](#); [Manual](#))

Clearing adapted from DOI: 10.1126/sciadv.aba0365

## Table of Contents

- [Table of Contents](#)
- [Setup](#)
  - [Prepare Stock Solutions](#)
    - [PBS-T \(1%\)](#)
    - [4% BSA in PBS](#)
    - [F-ara-EdU \(20mM stock in DMSO\)](#)
    - [Alexa Fluor 647 azide \(Component B\)](#)
    - [Click-iT reaction buffer \(Component D\)](#)
    - [Click-iT buffer additive \(Component F\)](#)
  - [Materials](#)
  - [Buffers](#)
- [Protocol](#)
  - [1. F-ara-EdU Pulse Labeling](#)
  - [2. Fixation & Permeabilization](#)
  - [3. Blocking](#)
  - [4. Click-iT F-ara-EdU Fluorescent Labeling and Detection](#)
    - [For 250µl reactions:](#)
  - [5. Staining DNA](#)
  - [6. Mounting and Imaging](#)
- [Notes/Ideas for Future Tweaking:](#)
- [Experiments](#)
  - [F-ara-EdU short chase](#)
  - [± chorion](#)
    - [Question: Does the chorion affect EdU incorporation?](#)

## Setup

### Prepare Stock Solutions

#### PBS-T (1%)

Reagent	Quantity (for 500 mL)	Final concentration
PBS 1X	495mL	-
Tween-20	5mL	1%

#### 4% BSA in PBS

- Dissolve 2g of BSA in 50mL of PBS.

## F-ara-EdU (20mM stock in DMSO)

- F-ara-EdU (2'S)-2'-Deoxy-2'-fluoro-5-ethynyluridine, (Sigma #T511293; FW 270.22 g/mol)
  - *Note: technically, given the formula weight, this is an 18.5mM stock...*
  - Add 1 ml of DMSO or PBS to vial containing 5 mg F-ara-EdU to make a 20 mM stock solution.
  - After being dissolved, remaining F-ara-EdU stock solution is stable for up to 1 year when stored at  $-20^{\circ}\text{C}$ .
  - F-ara-EdU solution should be portioned into single use aliquots to avoid freeze thawing.
  - CAUTION: EdU is incorporated into DNA and is a potential mutagen. Proper protective clothing should be used when handling EdU. Waste, including stock solutions, used substrate, and water containing F-ara-EdU should be considered as hazardous.

## Click-It stocks

### Alexa Fluor 647 azide (Component B)

1. Prepare a master stock solution of the Alexa Fluor® azide (Component B): Add 70  $\mu\text{L}$  of DMSO (Component C) to Component B, then mix well.
- After use, store any remaining working solution at  $\leq -20^{\circ}\text{C}$ . When stored as directed, this working solution is stable for up to 1 year.
2. Prepare a *Dye Working Stock* by diluting the master stock 1:100 in DMSO (assuming that the final concentration for the master stock from the kit is 5uM)

### Click-iT reaction buffer (Component D)

- Prepare a working solution of 1X Click-iT® EdU reaction buffer (Component D):
  - Transfer the solution (4 mL) in the Component D bottle to 36 mL of deionized water.
  - To make smaller amounts of 1X Click-iT® EdU reaction buffer, dilute volumes from the Component D bottle 1:10 with deionized water. After use, store any remaining 1X solution at  $2-8^{\circ}\text{C}$ . When stored as directed, this 1X solution is stable for 6 months.

### Click-iT buffer additive (Component F)

- To make a 10X stock solution of the Click-iT® EdU buffer additive (Component F): Add
  - 2 mL deionized water to the vial, then mix until fully dissolved.
- After use, store any remaining stock solution at  $\leq -20^{\circ}\text{C}$ .
- When stored as directed, this stock solution is stable for up to 1 year.
- If the solution develops a brown color, it has degraded and should be discarded.

## Custom Click stocks

### Azide Dye Stock

- picolyl-Azide-Dye: 4 mM in H<sub>2</sub>O (3 mg in 1 ml H<sub>2</sub>O) (final 8  $\mu\text{M}$ )

### Copper Sulfate Stock (CuSO<sub>4</sub>·5H<sub>2</sub>O; MW 249.68)

- Dissolve 250 mg  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in 5 ml  $\text{H}_2\text{O}$  for 200mM working stock
- Final concentration in reaction: 2mM

### Sodium Ascorbate Working Stock (Buffer additive, $\text{NaC}_6\text{H}_7\text{O}_6$ ; MW 198.11)

- **NOTE: Make Fresh for each experiment**
- Dissolve 0.198 g in 1 mL  $\text{H}_2\text{O}$  for 100mM working stock
- Final concentration in reaction: 10 mM

### Clearing stocks

#### Detergent Mix (1 L)

Reagent	Quantity (for 500 mL)	Final concentration
10% SDS	100mL	1%

- 1% SDS: 100 mL 10% SDS (may be filtered)
- 0.5% Deoxycholate: 5g Deoxycholate (wear a mask when handling the powder) 50 mM Tris-HCl (pH 7.5): 50 mL 1 M Tris-HCl (pH 8.0) 1 mM EDTA (pH 8.0): 2 mL 0.5 M EDTA (pH 8.0) 150 mM NaCl: 30 mL 5 M NaCl

#### Clearing Solution-1.1

- 8 to 10% (v/v) THEED (Sigma-Aldrich, 87600-100ML)
- 5% (v/v) Triton X-100 (Roth, 3051.2), and
- 5% (w/v) urea (Roth, X999.2) in  $\text{dH}_2\text{O}$

### Materials

- 1.5ml Microcentrifuge tubes (clear)
- Rocker
- Staining vials
- Click-iT EdU Alexa Fluor 647 Imaging Kit (ThermoFisher #C10340)

## Protocol

### 1. F-ara-EdU Pulse Labeling

1. Collect embryos at the desired stage.
2. Dechorionate embryos in Filtered Natural Sea Water or ASW.
  - TODO: Test whether dechoriation is necessary. If not, will permit more pulse and longer chases throughout development.
3. Transfer embryos or larvae into 1.5ml microcentrifuge tube or multi-well plate (e.g. 24-well plate) depending on desired pulse length.
  - If using a plate, you should pre-coat with 0.2% agarose to prevent embryos from sticking and allow for chorion expansion.

4. To optimize labeling efficiency and minimize toxicity, you should set up a test matrix of different concentrations and durations of F-ara-EdU labeling.

- Include cultures treated with vehicle alone for controls.

	30 min	60 min	120 min	24 h
1μM				
10μM				
100μM				
0μM (no EdU control)	x	x	x	x

5. Add appropriate volume of F-ara-EdU (20mM stock) to achieve desired final concentration.

- For example, to achieve a final concentration of 10μM, add 1μl of 20mM F-ara-EdU stock solution to 2ml of FNSW.

6. Incubate embryos in F-ara-EdU solution at 18°C for desired time interval.

7. If chasing, change back into FNSF + 1% Pen/Strep after pulse for duration of chase interval.

- Transfer to a clean plate or tube to avoid carryover of f-ara-EdU (e.g. in agarose).

8. Proceed immediately to fixation and permeabilization (below).

2. Fixation & Permeabilization

1. Fix embryos overnight in 4% PFA in FNSW at 4°C in round bottom 2mL tubes.
  - Alternatively, you **may** fix embryos for ~1h at RT in 4% PFA in FNSW but O/N is better.
  - Can be left for up to 4 days?
2. After fixation, remove fixation buffer and wash embryos in 1mL PBS-T (1%) on rocker.
  - 1x 5 min
  - 2x 10 min
3. Proceed to next step.

3. Clearing (optional)

1. Remove PBS-T and add 1mL of prechilled acetone (-20°C) to each tube.
2. Keep at -20°C for a minimum of 2h.
3. Wash 3x with 1mL PBS-T (1%) for 5-10 min each.
4. OPTIONAL (not recommended in *conjunction* with detergent below) Add 500μl of proteinase K (75 μg/mL; Merck 1245680100) and rock for 5 min at RT.
  - Quench with 2x exchanges of 1mL glycine (2mg/mL)
5. **PREFERRED Wash 3x with 1mL Detergent mix for 20-30 min each at RT.**
6. Wash 3x with 1mL PBS-T (1%) for 5-10 min each.
7. Transfer embryos to glass vial containing 10mL of Clearing Solution-1.1
  - Gently rock for 3-6 hours at 37°C
8. Transfer back to 2mL round bottom tubes and wash 3x with 1mL PBS-T (1%) for 5-10 min each.

4. Blocking

- 1. Block embryos in 1mL 4% BSA in PBS for 1h at RT on rocker.
- 2. Wash embryos in 1mL 4% BSA in PBS 1 x 5 min on rocker.
- 3. Proceed to next step.

PAUSE POINT: Store at 4°C

5. Click-iT F-ara-EdU Fluorescent Labeling and Detection

- 1. Allow Click-iT kit components to come to room temperature before opening.
  - 2. Prepare **1X Click-iT® EdU buffer additive** (Component F) by diluting the 10X stock solution 1:10 in deionized water. Prepare this solution fresh and use the solution on the same day.
  - 3. Prepare **Click-iT® reaction cocktail** according to Table. It is important to add the ingredients in the order listed in the table; otherwise, the reaction will not proceed optimally. Use the Click-iT® reaction cocktail within 15 minutes of preparation.
- For each tube mix together:
    - **430 µl of 1x Click-iT reaction buffer**
    - **20 µl CuSO4** (Component E)
    - **1.2 µl fluorescent azide** (Component B)
    - **50 µl 1x reaction buffer additive** (Component F; diluted above from 10x stock).
  - Scale up the mixture for the number of samples to be treated and add 500 µl to each tube.
  - It is important to use the cocktail within 15 min of preparation.
  - It is good practice to include a control sample of cells not exposed to EdU. In addition, these cells are needed for single staining compensation controls for intracellular antigens or antigens stained with RPE, PE-tandem, or Qdot antibody conjugates.

Table 3. Click-iT® reaction cocktails.

Reaction components*	Number of coverslips						
	1	2	4	5	10	25	50
1X Click-iT® reaction buffer (prepared in step 1.4)	430 µL	860 µL	1.8 mL	2.2 mL	4.3 mL	10.7 mL	21.4 mL
CuSO <sub>4</sub> (Component E)	20 µL	40 µL	80 µL	100 µL	200 µL	500 µL	1 mL
Alexa Fluor® azide (prepared in step 1.3)	1.2 µL	2.5 µL	5 µL	6 µL	12.5 µL	31 µL	62 µL
Reaction buffer additive (prepared in step 4.1)	50 µL	100 µL	200 µL	250 µL	500 µL	1.25 mL	2.5 mL
Total volume	500 µL	1 mL	2 mL	2.5 mL	5 mL	12.5 mL	25 mL
<b>*Note:</b> Add the ingredients in the order listed in the table.							

5. ALTERNATE Custom Click Rxn Table

Reaction Components / Reaction	1	2	4	5	10	25	50
1x PBS	394 µl	x µl	x µl	x ml	x ml	x ml	x ml
CuSO4 (200 mM)	5 µl	x µl	x µl	x µl	x µl	x µl	x µl
Fluorescent azide (4 mM)	1 µl	x µl	x µl	x µl	x µl	x µl	x µl

Reaction Components / Reaction	1	2	4	5	10	25	50
DMSO	50 $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x ml
Sodium Ascorbate (100 mM; Fresh)	50 $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x $\mu$ l	x ml
Total volume	500 $\mu$ l	x $\mu$ l	x ml	x ml	x ml	x ml	x ml

\* Add reagents in the order listed above. Use within 15 min of preparation.

\* DMSO protects the DNA from free radicals generated by the Cu(I) catalyst.

\* For lower dye concentrations, you can try 10  $\mu$ l of 10 $\mu$ M dye in 500  $\mu$ l reaction for a final conc. of 200nM.

1. Remove the blocking buffer.
2. Quick wash each sample *twice* with 1 mL of PBS. Remove the wash solution.
3. Add 0.5 mL of Click reaction cocktail to each tube. place on rocker.
4. Incubate for 30-60 minutes at room temperature, protected from light (longer gets deeper in hatchlings).
5. Remove the reaction cocktail and wash twice with 1 mL of 4% BSA in PBS.
6. Optional: Keep overnight in 4% BSA in PBS at 4°C to allow excess dye to diffuse out of yolk.

## 6. Staining DNA

1. Add 1 $\mu$ l of 1000X DAPI (**TODO: Need concentration**) to each tube.
2. Stain for 30-60 min at RT on rocker. (longer gets deeper into hatchlings).
3. Remove DAPI and wash once with 1 mL of PBS w/ 4% BSA.
4. Final wash in 1x PBS
5. Store at 4°C until imaging.

## 7. Mounting and Imaging

1. Mount on Matex in 1% agarose in Fluoromount-G
2. Image on confocal microscope.

## Notes/Ideas for Future Tweaking:

- Background in yolk significantly reduced with pre-block with BSA.
  - May also try casein (milk) for blocking too. or LMW-Dextran.
- A few suggestions for optimizing:
- Can we reduce the amount of fluorescent dye?
- Longer washes after Click (trying ON soak in PBS currently)
- Wash in high DMSO ( $\leq 50\%$  w/v) to remove unbound dye?
- 5% Urea wash to disrupt hydrophobic interactions with yolk?
- Can likely scale down the Click Rxn to 250 $\mu$ l in 2mL tubes. (This did not work so well first time.)
- Probably need to consider:
  - Increasing the dose  $>100\mu$ M for shorter pulses
    - Need to read up on what is usual dosing for EdU in other systems (e.g. zebrafish embryos; 400  $\mu$ M EdU for 1 h at room temperature?)

- Increasing the pulse time (1 hour is not really labeling much, and 24 hours seems good, but maybe because some cells are cycling twice?)
- 

## Experiments

---

### F-ara-EdU short chase

- Pulse St20 embryos for 1 hr with 100μM F-ara-Edu
  - 5mL FNSW required per dish
    - 25μl 20mM F-ara-EdU per 5mL FNSW
- change back into FNSF + 1% Pen/Strep after pulse
  - Reserve some 0hr chase embryos for comparison
- Collect after:
  - 0h
  - 1h
  - 2h
  - 4h
  - 8h
  - 12h
  - 24h

### Results

- Need to wash embryos more thoroughly after Click Rxn
- May need higher dose (400μM?) and/or longer pulse for better labeling

± chorion

### Question: Does the chorion affect EdU incorporation?

- Pulsed St26 embryos ± chorion in 50μM F-ara-Edu (12.5μl per 5ml FNSF) for:
  - 1h
  - 24h
- Blocked for 1h
- Only used 250μl Click Rxn
- Result for 1 hr is very weak/low staining with moderately high background:
  - likely because of less proliferation at later stage, short pulse, and/or accessibility
- Need to assess the 24 hr pulse as well
  - 24 hour pulse click Rxn performed in 500 mL reaction
  - 24 hour pulse much stronger and cleaner

### Results

- Lower levels of F-ara-EdU incorporation in embryos with chorion
- Interesting structures observed on anterior side.

## Suggestions from the Lab

- Bring dye concentration down to ~50nM (vs 5 $\mu$ M final in standard ClickIt kit)
  - Azdye 647 - azide plus (Click chemistry tools)
  - 50nM
  - On ice, one minute
  - Get Richard's protocol for EdU-click
    - 1 mM CuSO<sub>4</sub>
    - 200 nM AzPlus-dye in PBS-T w/ 10% DMSO (This is 2 orders of magnitude lower than the Alexa647 in the ClickIt kit)
    - 5 mM sodium ascorbate
  - Stocks:
    - 100mM CuSO<sub>4</sub>
    - 100mM sodium Ascorbate